



MDS1  
Middle School Microscope



## **Magnification Options**

Supplied with a Widefield 10x eyepiece, along with 4x, 10x, 40x (retractable) objectives to provide a range of magnification options: 40x, 100x, 400x.

## **Ease of Use**

Coarse and Fine focus knobs make it easy to obtain a crisp and clearly focused image. Rack stop protects objectives and specimen slides from damage. Head rotates 360°.

Specimen circular stage slides forward, backward and left-right. Locked-on spring-activated stage clips hold samples in place.

Rugged, all-metal frame has enamel coating for protection. Carrying handle protects microscope from damage during transport.

## **Cordless or Corded LED Illumination**

Cool, bright LED Illumination. Includes rechargeable batteries and recharger. Use corded or cordless. CE certified power cord. 100~240V. Rheostat control for illumination.

## **You're Covered!**

Includes dust cover and 5-year warranty.



MDS1 Microscope Magnification

Achromat Objective	10x Eyepiece
4x	40x
10x	100x
40x	400x