

S2 Microscopes





S2 Specifications

- Dual Magnification
- Objective Lenses: 1x & 3x
- Eyepieces: WF10x/20, WF15x/15, WF20x/10
- Diopter adjustment on left eyetube
- Working Distance: 100mm
- 3 year limited warranty



Model: S2-SPS-LCD



Model: S2-SPS with Ring Light

S2 Models

S2-DL

Binocular or Trinocular, Track Stand with top and bottom light.

S2-SPS

Binocular or Trinocular, Post Stand

S2D-SPS

Trinocular, C-Mount Adapter, 5.1mp Digital Camera, Post Stand, LED Ring Light.

S2-SPS-LCD

Trinocular, C-Mount Adapter, 5mp Android Tablet Camera, Post Stand with LED Ring Light.



Model: S2-DL

info@richter-optica.com
www.richter-optica.com