



### Color Camera Specifications

- USB-C + WiFi Output
- 8 Megapixels
- 4K HD Camera
- Sensor Size: 1/1.8" back-illuminated CMOS Sony IMX678 Starvis 2 sensor
- 30 frames per second at 3840 x 2160
- Pixel Size: 2.0 x 2.0 $\mu$ m
- Rolling Shutter
- Capture images with PC Software or App
- PC Software allows for measurement and annotation
- 3D Noise Reduction
- ROI White Balance
- Anti-Flicker
- C-Mount connection to microscope
- Size & Weight: 1.97" x 1.97" x 4.4", 1.1 lb.
- Included: Camera, Antenna, USB cable, US plug.

### Applications:

- Education
- 4K Microscopy Imaging
- Routine Industrial Inspection
- Exhibit Display
- Remote Observation

