

## S - 2

### Stereo Microscope

#### Specifications

- Dual Magnification
- Objective Lenses: 1x & 3x
- Eyepieces: WF10x/20x, WF15x/15, WF20x/10
- Diopter adjustment of left eyetube
- Working Distance: 100mm
- 3 limited warranty



Model: S2-SPS-LCD



Model: S2-SPS with Ring Light

#### Models

##### S2-BL

*Binocular or Trinocular, Track Stand with top and bottom light*

##### S2-SPS

*Binocular or Trinocular, Post Stand*

##### S2D-SPS

*Trinocular, C-Mount Adapter, 3mp Digital Camera, Post Stand, LED Ring Light*

##### S2-SPS-LCD

*Trinocular, C-Mount Adapter, 5mp Android Tablet Camera, Post Stand with LED ring light*



Model: S2-BL